

Issue Date: 2009-04-10  
2010-07-03

Page 1 of 9

Report Reference #

E248249-A4-UL

## UL TEST REPORT AND PROCEDURE

<b>Standard:</b>	UL 61010-1, 2nd Edition, 2005-07-22 (Electrical Equipment for Measurement, Control, and Laboratory Use; Part 1: General Requirements) CAN/CSA-C22.2 No. 61010-1, 2nd Edition, 2004-07, (Electrical Equipment for Measurement, Control, and Laboratory Use; Part 1: General Requirements)
<b>Certification Type:</b>	Laboratory Use Electrical Equipment
<b>CCN:</b>	OGTK, OGTK7
<b>Product:</b>	AIRSTREAM VERTICAL LAMINAR FLOW CLEAN BENCH
<b>Model:</b>	AVC-2D2, AVC-3D2, AVC-4D2, AVC-5D2 and AVC-6D2
<b>Rating:</b>	110 - 130 Vac, 60 Hz AVC-2D2 350W/3.5A AVC-3D2 380W/3.5A AVC-4D2 400W/3.5A AVC-5D2 450W/3.5A AVC-6D2 800W/6.5A
<b>Applicant Name and Address:</b>	ESCO MICRO PTE LTD 21 CHANGI SOUTH ST 1 SINGAPORE 486777 SINGAPORE

This is to certify that representative samples of the products covered by this Test Report have been investigated in accordance with the above referenced Standards. The products have been found to comply with the requirements covering the category and the products are judged to be eligible for Follow-Up Service under the indicated Test Procedure. The manufacturer is authorized to use the UL Mark on such products which comply with this Test Report and any other applicable requirements of Underwriters Laboratories Inc. ('UL') in accordance with the Follow-Up Service Agreement. Only those products which properly bear the UL Mark are considered as being covered by UL's Follow-Up Service under the indicated Test Procedure.

The applicant is authorized to reproduce the referenced Test Report provided it is reproduced in its entirety.

Any information and documentation involving UL Mark services are provided on behalf of Underwriters Laboratories Inc. (UL) or any authorized licensee of UL.

Prepared by: Lesley Tee  
Underwriters Laboratories Inc.



Reviewed by: Chanyo Won  
Underwriters Laboratories Inc.



Copyright © 2010